

A-dec® Choices

Which Solution is Right For You?



a dec®

There's an A-dec for everyone.

A-dec has collaborated with the dental community for decades, ensuring that our equipment is designed to work the way you do. And while you may not keep a 20-year-old chair in your operatory, rest assured, it's not uncommon for A-dec equipment to transfer from generation to generation. No matter which chair you choose, you'll have peace of mind knowing that we build our equipment like our reputation depended on it. Because it does.



A-dec 500 Optimal access. Flexible integration. Intelligent control. The A-dec 500 brings intelligent simplicity and infallible peace of mind to your operatory, and superior comfort to your patients.



A-dec 400 Sleek design, patient comfort, and uninhibited dental team access combined with solid, reliable performance.



A-dec 300 Stylish and compact. A-dec 300 combines a robust design, great access, and minimal maintenance.



A-dec 200 Part of a complete system packed with features for added accessibility and comfort—all at a great value.

Chair	A-dec 500	A-dec 400	A-dec 300	A-dec 200
Ultra-soft start/stop hydraulic drive	•			(without soft start/stop)
Hybrid drive (hydraulic base/electric tilt)		•	•	
Swivel (30° each side of center)	•	•		
Virtual pivot for added patient comfort	•	•	•	
Thin backrest 32mm	•	•	•	
Dual-articulating gliding headrest	Lever release	Locking knob or Lever release	Locking knob	Locking knob (without glide)
Multi-position armrest settings	2	2	2	2
Chair at lowest point	343mm	349mm	349mm	381mm
Chair at highest point	800mm*	800mm	800mm	737mm
Integrated 300 watt power supply	•	•	•	•
Load rating	227 kg.	181 kg.	181 kg.	181 kg.
Choice of sewn or formed upholstery	•	•	•	Formed only
Delivery System				
Delivery system	A-dec 500	A-dec 300	A-dec 300	A-dec 300, A-dec 200
Four-position control block	Standard	Optional	Optional	Standard
Three-position control block	—	Standard	Standard	—
Holder positions	6	5	5	5
Clinical device integration	5	2	2	Yes
Deluxe plus touchpad	•			
Choice of programmable standard or deluxe touchpad		•	•	•
Support center	A-dec 400	A-dec 400, A-dec 300	A-dec 400, A-dec 300	A-dec 200
Self-contained water system with quick-disconnect 2-liter bottle	•	•	•	(without quick-disconnect)
Accommodates amalgam separator	•	•	•	•
Cuspidor swivel	•	(not on the fixed glass cuspidor)	(not on the fixed glass cuspidor)	•
Programmable timing for cup and bowl	•	•	•	•
Dental light	A-dec 500 LED	A-dec 500 LED, A-dec 300 LED	A-dec 500 LED, A-dec 300 LED	A-dec 200, A-dec 500 LED, A-dec 300 LED

*Using maximum upper limit. The factory set upper limit is 737mm.



A-dec Worldwide

With A-dec, you gain a partner for the life of your practice. From our headquarters in Newberg, Oregon, we support global offices and a network of dedicated dealers in more than 100 countries worldwide.

Discover the A-dec Difference. Visit a-dec.com/TheDifference to find out how every detail behind an A-dec solution ensures a reliable, lasting investment.



A-dec Headquarters

2601 Crestview Drive
Newberg, Oregon 97132 USA
Tel: 1.800.547.1883 within USA/CAN
Tel: +1.503.538.7478 outside USA/CAN
Fax: +1.503.538.0276
a-dec.com

A-dec Australia

Tel: 1.800.225.010 within AUS
Tel: +61.(0)2.8332.4000 outside AUS

A-dec China

Tel: 400.600.5434 within CN
Tel: +86.571.89026088 outside CN

A-dec United Kingdom

Tel: 0800.233285 within UK
Tel: +44.(0).24.7635.0901 outside UK



Certain touchpad symbols are proprietary to A-dec Inc. Any use of these symbols, in whole or in part, without the express written consent of A-dec Inc. is strictly prohibited.

The colors and products shown may vary from the actual materials and are subject to change without notice. Please see your authorized A-dec dealer for accurate color samples and the most current product information.

©2022 A-dec Inc. All rights reserved.
85.0394.25/E-only/11-21(2-22)/Rev C2