

# Monitor Mounts



# COMPLETE PICTURE. CLEAR VIEW. DEPENDABLE OPTIONS.

No matter which option you choose, you and your team will enjoy an ergonomic solution that offers richer patient-staff interactions. Robustly built and painstakingly tested, all monitor mounts are VESA compatible

to accommodate various sizes\* of medical-grade monitors up to 20 pounds (9 kg).

## CHAIRSIDE: INCREASE CASE ACCEPTANCE.

A chairside monitor brings the screen closer to the patient, so it's easier to discuss oral health, explain treatment procedures, and gain case acceptance.

CHAIRSIDE MONITOR MOUNT



- **Close-to-the-patient positioning.** Chairside mount encourages consultation and education.
- **Tilt and swivel.** Offers comfortable viewing even in seated and supine positions.
- **Easy-grip handle.** Maneuver and tilt the monitor with one hand.

## INTEGRATED: GET A BETTER VIEW.

Integrate a monitor into A-dec Inspire® cabinets. The A-dec gliding monitor mount slides and pivots for picture-perfect positioning between you and your assistant. Mounted at 12 o'clock, it's out of the patient's direct view.

A-DEC INSPIRE® GLIDING MONITOR MOUNT



- **Smooth and precise.** Telescoping monitor glides vertically and horizontally for task-based positioning.
- **Clean and sleek.** Design neatly conceals cables inside the arm, out of site.
- **Positionable for privacy.** At 12 o'clock, the monitor can move outside of treatment room for discussions

\* Ensure monitor's dimensional size does not interfere with other equipment.

## CEILING/WALL/CABINET: EXPAND VISIBILITY.

Keep the monitor within reach when you need it, and move out of the way when you don't.

- **Exceptional range of motion.** Get a full spectrum of viewing angles for you, your team, and your patients.
- **Versatile flexarm design.** Easy adjustments and smart cable routing simplifies installation and streamlines appearance.
- **Space-saving design.** Moves out of the way during procedures, and works well for patient entertainment.

CEILING MOUNT




CENTRAL CABINET MOUNT



WALL/SIDE CABINET MOUNT



# WHAT POSITION WORKS BEST?



## A-dec Headquarters

2601 Crestview Drive  
Newberg, Oregon 97132 USA  
T: 1.800.547.1883 within USA/CAN  
T: +1.503.538.7478 outside USA/CAN  
F: +1.503.538.0276  
a-dec.com

## A-dec Australia

T: 1.800.225.010 within AUS  
T: +61.(0)2.8332.4000 outside AUS

## A-dec China

T: 400.600.5434 within CN  
T: +86.571.89026088 outside CN

## A-dec United Kingdom

T: 0800.233285 within UK  
T: +44.(0).24.7635.0901 outside UK

All trademarks in this document are subject to the trademark rights of A-dec, its affiliates, and their collective suppliers in the United States and foreign countries.

The colors and products shown may vary from the actual materials and are subject to change without notice. Please see your authorized A-dec dealer for accurate color samples and the most current product information.