

# Integration has never been so easy.

---

Now, A-dec® delivery systems with 300 Deluxe and 500 Deluxe Plus touchpads accommodate the KaVo® KL 703 LED electric micromotor. The streamlined design and intuitive user interface make the system easy to operate.

Enjoy an exceptional practice experience, with premier equipment that complements your professional expertise and reflects your commitment to superior patient care. Elevate your practice with A-dec and KaVo.



**KaVo KL 703 Electric Motor: Smart. Efficient. Streamlined.**



# INTRA™ LUX KL 703 LED Electric Motor



## LIGHTWEIGHT FEEL. COMPACT DESIGN. ERGONOMIC SHAPE.

Comfort and performance come together in the KaVo KL 703. Lightweight and compact, this sterilizable motor provides constant torque over the entire speed range for fatigue-free comfort and endless power. Suitable for any A-dec 200, A-dec 300 and A-dec 500 delivery system with a Deluxe or Deluxe Plus touchpad.

### Features and Benefits

- Seamless electric motor integration with the choice between single and dual motors
- High-speed, low-speed, and torque-controlled rotary and reciprocating endo applications
- Speed range between 100 and 40,000 rpm with a 1:1 transmission ratio
- Brushless and sterilizable
- Up to 8 speed memory presents available in prep mode
- Up to 16 programmable memory settings (with A-dec 500 Deluxe Plus touchpad)
- ISO 6964 E-style coupling

<b>Applications</b>	General restorative, prophylaxis, endodontics	<b>Length (mm)</b>	34
<b>A-dec Delivery System</b>	A-dec 500B / A-dec 300 / A-dec 200	<b>Weight (gm)</b>	69
<b>System</b>	E-style	<b>Sterilizable</b>	Yes
<b>Power (watts)</b>	45	<b>Micromotor/Type</b>	Brushless motor
<b>Speed (rpm)</b>	100 to 40,000	<b>Warranty</b>	3 years
<b>Torque (Ncm)</b>	Max 3.0		

## ELECTRIC MOTOR SYSTEM

Control Module + Motor + Tubing Kit			Motor + Tubing Kit	
<b>Control Module</b>			<b>Motor</b>	
	EMC-703 electric motor controller			
<b>Motor</b>			<b>Tubing Kit</b>	
	INTRA LUX KL 703 LED electric motor			
			Electric motor tubing (173 cm medium or 208 cm long)	
<b>1st Position</b>			<b>2nd Position</b>	